

Corrosiveness to Copper from Petroleum Products

ASTM D130 - EN/ISO 2160 - IP 154 - DIN 51811 - ASTM D4048 - IP 112 - ASTM D7095



Stainless steel bath
Cooling coil
Bath drain
Easy to operate
Position for 6 test cylinders (9 optional)

Item	Unit	TC16
Ordering code		
TC16 230V Circulator	50 or 60Hz	31T0671
TC16 115V Circulator	60Hz	31T0861
Power	[kW]	1.5
Range	°C °F	Ambient .. 250 Ambient ..482
Reading		°C or °F
Setting	[°]	0.1
Stability	[°K]	0.02
Heating	[kW]	1.4
Bath volume	[L]	16
Bath opening	[mm]	6 x 51 (9 x 51 optional)
Bath depth	[mm]	220
Length	[mm]	480
Width	[mm]	295
Height	[mm]	480
Materials	Used inside bath: stainless steel 304, brass	
CE	Conforms to CE regulation	

General

The test method covers the detection of the corrosiveness to copper of aviation gasoline, aviation turbine fuel, automotive gasoline, natural gasoline or other hydrocarbons having a Reid Vapour Pressure no greater than 18 psi. The cover of the bath has six openings, each with a lid and hook for suspending a copper corrosion test vessel (P/N 14T0100). Or the openings can accommodate test tubes in holders (P/N 14T0102). The temperature range is from ambient +5°C to 250°C. Optional is a cover with nine holes to test nine samples (P/N 03T2311).

Accuracy

The insulation of the bath and electronic design result in a very stable working temperature of $\pm 0.02^\circ$. The set point can be set in steps of 0.1° in the range of 0°C up to 250°C ($-148..482^\circ\text{F}$). The accuracy on the display is displayed in 0.1°C . However the controller has an internal accuracy of 0.01°C .

Temperature readout

Standard available in $^\circ\text{C}$, on request in $^\circ\text{F}$.

Pump

When not used for copper corrosion tests, the pump can be used to circulate the bath content to an external application.

Safety

The bath conforms to CE regulation. It is further equipped with a mechanical resettable safety thermostat. Optional is a low liquid floater (P/N 07T0080).










Alternative

We also offer a TC40 bath with a cover for 18 positions (P/N 03T2313).









Corrosiveness to Copper from Petroleum Products

Accessories

P/N	Picture	Description
14T0100		Pressure vessel for copper strip corrosion test, test cylinder with o-ring. Cylinder is supplied with works certificate.
14T0101		Lid and mounting hook (six pieces standard included)
13T8000		Top lid with six 51ø mm holes and thermometer hole (standard included)
31T0003		ASTM copper strip corrosion test standard
31T0002		Copper strip (per one strip) 12.5 x 3.0 x 75 mm (width x thick x length) 99.9 % purity
31T0001		Sanding paper silicon carbide P220
09T0010		Test tube 25 x150mm







Corrosiveness to Copper from Petroleum Products

Accessories

P/N	Picture	Description
14T0102		Test tube (25 x 150 mm) holder
31T0000		Multistrip vise, holds up to four strips while polishing
03T2311		Top lid with nine 51ø mm holes. Please separately order 3* P/N 14T0101 additionally.
31T0005		Silicon carbide powder 105 µm (1 kg)
09T0011		Flat viewing test tube
08T0001		Bath fluid silicon oil 20..150°C 20ltr
25T0911		ASTM 12C thermometer
25T2154		Thermometer holder, 425 x 10mm
31T0008		Vented stopper made from silicon rubber for test tube (P/N 09T0010). 25 pieces. Width 24 (top) to 18 (bottom) mm. Length 30 mm. Hole of 4 mm. Temperature resistant from -60 to +180°C.
31T0009		Sanding paper silicon carbide P240

Corrosiveness to copper from Petroleum Products

Accessories




P/N	Picture	Description
24T0385		O-ring for test cylinder (P/N 14T0100)
31T0010		Disposable copper foil strip for ASTM D7095 (one piece)
31T0011		Scouring pad 400 grit
14T0109		Test tube holder test jar conform ASTM D4048
31T0012		Test jar ASTM D4048
31T0013		Vented cork for test jar (P/N 31T0012). 10 pieces.

Corrosiveness to copper from Petroleum Products

Accessories

P/N	Description	ASTM D130, IP 154	IP 112, ASTM D4048	ASTM D7095
14T0100	Test vessel with o-ring	○		○
31T0003	Test standard	○	○	○
31T0002	Copper strip	○	○	
31T0001	Sanding paper P220	○		
09T0010	Test tube 25 x150mm	○		○
14T0102	Test tube holder	○		○
31T0005	Silicon carbide powder	○	○	
31T0000	Multistrip vise	○	○	○
09T0011	Flat glass viewing testtube	○	○	○
25T0911	ASTM 12C thermometer	○	○	○
08T0001	Bath fluid silicon oil 20-150°C 20ltr		○	
25T2154	Thermometer holder	○	○	○
31T0009	Sanding paper P240		○	
31T0008	Vented cork for test tube	○	○	○
24T0385	O-ring for test vessel (P/N 14T0100)	○		○
31T0011	Scouring pad 400 grit			○
31T0010	Disposable copper foil strip			○
14T0109	Test jar holder ASTM D4048		○	
31T0012	Test jar ASTM D4048		○	
31T0013	Vented cork for test for ASTM D4048		○	

Alternative setup offering 18 positions (instead of 9)

P/N	Picture	Description
00T0681		TC40 230V/ 50-60Hz
00T0851		TC40 115V/60Hz
03T2313		Lid with 18* 51ø mm holes and lids
14T0101		Lid + mounting hook number of pieces to be ordered separately depending on number of test cylinder used by customer